



While you focus on patients,
we're focusing on **you.**

Surgical Solutions

Always by your side

With a unique combination of infection protection and contamination control products and services, STERIS offers customized solutions for a variety of environments. Rely on STERIS to help you gain productivity, lower costs, and protect staff and patients throughout your OR, sterile processing department, GI/endoscopy suite, and ambulatory surgery center.

When it's a matter of life and death, patients and their families count on the expertise of caring professionals, who count on STERIS. Our mission is to help ensure that the environment in hospitals and surgical centers allows every patient to receive vital care without undue risk from infection or contamination. Our offering includes general and specialty surgical tables, surgical and examination lights, equipment management systems, and other complementary products and accessories. Our solutions give you the freedom to focus on what you do best – care for patients.



EXPERTISE &
EXPERIENCE




We've been serving you in the OR for decades.

We're dedicated to bringing you more of what you need at every point in and around the OR.

- Full range of solutions for decontamination and sterilization
- Versatile and user-friendly perioperative equipment
- Simplicity and sterility of single-use endoscopy devices
- Design, planning and installation of infrastructure for the OR and other procedural areas
- Innovative solutions for instrument management
- Comprehensive clinical education and consultation
- Knowledgeable technical service and clinical support
- Meaningful collaboration and proactive relationship management

We put our quality, reliability, service, and reputation to work for you every day.



Life in the OR is complex. We help you simplify it.

We understand your need to minimize complexity and focus on what matters most. Together, we help your team deliver uncompromised patient care.




Surgical tables and accessories

- General surgical tables: 4000, 5000 Series
- STERIS OT 1000 Series Orthopedic Surgical Table
- SurgiGraphic® 6000 Image Guided Surgical Table
- Surgical table accessories



Lighting

- HarmonyAIR® Surgical Lighting Systems
- HexaLux® Medical Examination Light
- Surgical lighting accessories



OR environment

- CLEANSUITE® OR Ceiling System
- MEDglas® Prefabricated OR Walls



OR integration systems

- HexaVue™
- HexaVue™ CONNECT
- Surgical displays and large format displays
- Utility and video connectivity systems




Equipment management

- HarmonyAIR® Equipment Column
- HarmonyAIR® Supply Head
- HarmonyAIR® Monitor Carrier
- HarmonyAIR® Smoke Evacuation System
- HarmonyAIR® Accessories



Case goods

- AMSCO® single and dual-compartment warming cabinets
- Scrub sinks
- Operating room storage solutions
- Nurse documentation stations



Meaningful, clinician-inspired advances

Keeping up with life in today's OR requires thoughtful innovation, bottom to top.



Support virtually any need

We offer a variety of reliable and versatile mobile, electro-hydraulically operated tables designed to support practically all surgical, MIS and robotic procedures, including cardiothoracic, endovascular, bariatric, robotic, ENT/ophthalmology, orthopedic, gynecologic/obstetric, and urologic/kidney.



Redefine optical performance

Our surgical lights feature excellent color, contrast, clarity, and control. They are designed for bright white light to illuminate the surgical site and eliminate shadows.



Improving your OR environment in measurable ways

We think about your environment with every solution we provide.



Purify

CLEANSUITE® OR Ceiling System air quality surpasses the standard for operating rooms by 1000X (US FED STD 209E CLASS 100,000/ISO8). It is also 6x faster to install and up to 40% less costly than traditional systems built onsite.

Strengthen

MEDglas® Prefabricated OR Walls combine beauty, strength, ultra-durability, excellent cleanability, and can withstand the force of a 9.2 magnitude earthquake.

Simplify

Our OR integration solutions allow for remote monitoring to minimize traffic in the OR. Displays are fanless and ventless for easier cleaning and a quieter atmosphere conducive to the work of healing.

Other thoughtful approaches such as booms, ConnectPoint, and Utility Hub keep cords and cables off the floor and out of your way for safety. Lighting features an airfoil arm design and cutout in the lighthead to promote laminar airflow. Even our patient warming system is conductive, not convective.

Every second matters



Integrate technologies seamlessly

We understand that the OR is a fast, dynamic environment that depends on efficiency, cleanliness, and safety. Our integration systems are built specifically for the OR and configured to meet your needs now and for years to come.

Control and manage the technology of an entire OR from one centralized access point. The issues that matter to you, such as cybersecurity and telemedicine, matter to us, too. The HexaVue platform was developed with an uncompromising focus on cybersecurity and a commitment to maintaining it.



The HexaVue integration system supports truly uncompromised 4K visualization and image capture.



Our integration products support you in the rapid shift to telemedicine.



Make every day easier

Storage that meets your needs helps every day run more smoothly. Warming cabinets meet AORN recommended practices for a safe environment of care, including a design that is purposefully built for warming IV solutions. Constructed entirely of stainless steel and designed in accordance with UL safety standards, our range of scrub sinks can fulfill the requirements of any operating room.

Storage solutions are designed for multiple specialties, and to adapt to new workflows. We can also provide an ergonomic and efficient workspace for reporting, charting, and management.



Maximize room space

Organize surgical equipment and medical services in one central location, delivering 360° of equipment accessibility. Our booms feature ultra-glide bearings to move smoothly and easily with little effort, while an electromechanical braking system offers reliable and drift-free movement.

Bringing it together from the start

We offer peace of mind and responsive, knowledgeable technical assistance and troubleshooting to support you in your mission of patient care. Our passionate pursuit of quality helps us develop technologies that deliver more to you and your patients, with flexible plans to meet any need and budget.



From the ground up

Our professionally staffed Project Design and Planning Group assists healthcare providers, architects, and facility planners in the design of hospital renovations and expansion. We help ensure that every space will be functional, practical, and thoughtfully designed to support healing.

This team of field-based project design managers, installation managers, and dedicated medical equipment planners has an average of 25 years of experience in hospital and healthcare facility planning and design. Over nearly six decades, the STERIS project design management group has completed thousands of projects worldwide and projects in more than half of the hospitals in the USA.



Always at your service

With hundreds of trained and experienced field service technicians in the USA, we provide service that's reliable, convenient, and only a phone call away. Your service provider not only has a backup, but also an additional back-up. Our bench of all-stars is waiting to serve you.



**Together, we're
building the OR
of your future.**



Future-ready to evolve with you

All of us here at STERIS are committed to you through our decades of experience, hundreds of service technicians in the USA, and an extensive line of products for hospitals and ambulatory surgery centers.

- **Surgical tables and accessories**
- **Lighting**
- **OR environment**
- **OR integration systems**
- **Equipment management**
- **Case goods**

For more information, contact your STERIS representative or visit www.steris.com

Document # M10824EN.2020-10, Rev.A
©2021 STERIS
All rights reserved.

 **STERIS**

STERIS
5960 Heisley Road
Mentor, OH 44060-1834 • USA
440-354-2600 • 800-548-4873
www.steris.com