

California Ingredient Disclosure



Kart-Klenz III Concentrated Liquid Detergent

Product Number: 1139

Revision Date: 01/01/2020

Section 1: General Information

Description of Product: Highly Concentrated Liquid Detergent

Manufacturer: STERIS Corporation

URL to SDS: <https://ww1.steris.com/onbDocs/V392/1213/645357.pdf>

Section 2: Functional Ingredients*

Chemical Name	CAS Registry Number	Functional Purpose	URL(s) to Designated List(s)
2-Propenoic acid, sodium salt, polymer with 2-methyl-2-((1-oxo-2-propenyl)amino)-1-propanesulfonic acid monosodium salt	37350-42-8	Corrosion and Scale Inhibitor	N/A
Etidronic Acid	2809-21-4	Chelating Agent Scale Inhibitor	N/A
Sodium Hydroxide	1310-73-2	Alkalinity Source pH Adjustment	N/A
Sodium Xylene Sulfonate	1300-72-7	Hydrotrope Coupler Surfactant	N/A
Sodium Sulfite	7757-83-7	Corrosion Inhibitor	N/A
Lauramine Oxide	1643-20-5	Surfactant	N/A
D&C Blue MX 176	N/A	Dye	https://monographs.iarc.fr/agents-classified-by-the-iarc/

California Ingredient Disclosure



Kart-Klenz III Concentrated Liquid Detergent

Product Number: 1139

Revision Date: 01/01/2020

Section 3: Non-Functional Ingredients*

Chemical Name	CAS Registry Number	URL(s) to Designated List(s)
N/A	N/A	• N/A

*Functional and Non-Functional Ingredients labeled at N/A do not meet listing requirements according to California Senate Bill Number 258.

Section 4: Additional Information

California SDS provided upon request by emailing: asksteris_msds@steris.com