


REVISION SCHEDULE		
REVISION #	DESCRIPTION	PLANNER/DATE



EQUIPMENT LIST - STANDARD EQUIPMENT

TAG	QTY	DESCRIPTION
01	1	STERIS Clean-Up Sink - 3 Well/120" Long
02	1	Caviwave Pro 17W Ultrasonic Cleaner - Wash/Rinse/17 Gallon
03	1	Acu-Hold System - Wall Mounted/2 Container
04	2	Ergotron LX Wall Mount System with 23" LCD Monitor
05	2	Reliance Vision Single-Chamber Washer/Disinfectant - Double Door/Steam Heat
06	4	Reliance Transfer Cart with Reliance 3 Level Manifold Rack
07	1	Amsco Pass-Thru Window with Set-Down Counter - Modular Wall
08	4	STERIS Ultra Frame Deluxe Single Prep and Pack Workstation - 30x60"
09	1	Amsco V-PRO Max Low Temperature Sterilization System - Single Door/Cabinet/Set of (4) Locking Casters

EQUIPMENT LIST - STANDARD EQUIPMENT (CONTINUED)

TAG	QTY	DESCRIPTION
10	2	Amsco 400 Medium Steam Sterilizer - 26x37.5x36" (660x953x910mm) (575L)/Double Left-Right Hinged Doors/Recessed 2 Walls/Steam Heat
11	4	Loading Car and Transfer Carriage - 36" (914mm)
12	2	Amsco Eagle Power Electric Steam Generator - 75 kW
13	1	AmeriWater Healthcare Water Treatment System
14	1	STERIS Blow-Down Separator Kit with Chromalox Tank
15	1	Emergency Exit Stairs and Platform (Required)
16	2	Hand Wash Sink With Eye Wash
17	1	Decontamination Room HVAC System
18	1	Sterilizer Cube HVAC System
19	1	Sterile Room HVAC System

EQUIPMENT LIST - OPTIONAL OR CUSTOMER SUPPLIED EQUIPMENT

TAG	QTY	DESCRIPTION
21	1	Decontamination Room Entrance Stairs and Platform (Optional Based on Deployment Configuration)
22	1	Sterile Room Exit Stairs and Platform (Optional Based on Deployment Configuration)
23	1	Amsco V-PRO Max Low Temperature Sterilization System - Single Door/Cabinet/Set of (4) Locking Casters - Optional (Customer Supplied)
24	2	6' Wire Storage Rack (Customer Supplied)