


Ultimate Savings

Smart materials management doesn't mean having to settle for less. **Prolystica® Cleaning Chemistries** allow you to use much less product in each cleaning cycle without sacrificing a single aspect of cleaning performance. Plus, Prolystica's patented formulations can increase instrument life and provide savings year after year.



“Instruments looked brand new that had been processed for two weeks.”¹

Save on instruments

Prolystica Cleaning Chemistries are specifically formulated to protect your instruments. Our patented formulas minimize the damaging effects of water as corrosion inhibitors prevent oxidation and pitting, resulting in longer instrument life.

Save on maintenance

Prolystica Cleaning Chemistries are also easier on your washers than the leading solid detergent system. The concentrated formulas can reduce maintenance costs and promote longer washer life.

Save on utilities

With Prolystica Cleaning Chemistries, you'll use less water and run fewer cleaning cycles for shorter periods of time. The low-foaming formulas eliminate the need for second rinses.

Save the planet

Prolystica Cleaning Chemistries are chemistries you can feel good about using, with proven biodegradability and recyclable packing, they're a smarter choice for your sterile processing department.

Instrument Protection Studies

Before



After only wash cycles




Compare after real application of Prolystica Ultra Concentrate Enzymatic Cleaner and Prolystica Ultra Concentrate Alkaline Detergent.

Prolystica Ultra Concentrate Enzymatic Cleaner 1/10 oz per gallon

Prolystica Ultra Concentrate Alkaline Detergent 1/10 oz per gallon

Scissor comparison after immersion in highest recommended use-dilution of each product in tap water for 5 weeks. First signs of rust noticeable in less than an hour on some products, and less than 24 hours on others.



Prolystica Ultra Concentrate HP
Neutral Detergent
1/10 oz per gallon

PowerCon®
Neutral pH
Detergent
1/2 oz per gallon

Tap Water

Tri-Power
Neutral
Detergent
1/10 oz per gallon

Asepti-Solid®
Neutral Detergent
1 oz per gallon
equivalent

SuperNova®
.1 NpHD
1/10 oz per gallon

¹Customer survey response at the end of independent STERIS product trial with Prolystica Ultra Concentrate HP Cleaning Chemistries.

Tri-Power is a product of United Biotech, SuperNova is a trademark of Case Medical, Inc., PowerCon is a trademark of Getinge USA, Inc., Asepti-Solid is a trademark of Ecclab, Inc.



Scan this code or visit www.steris.com/prolystica to learn more.

