




Integration for the Most Demanding Surgical Environments

Harmony® iQ 3600
OR Integration System



Supporting the Most Demanding, Complex ORs Now and Into the Future

STERIS recognizes the challenges facing today's ORs:

- Achieving efficient workflow through tools that are easy to use
- Keeping a focus on the patient to deliver excellent care instead of diagnosing OR equipment issues
- Connecting the devices, systems and teams needed to achieve better outcomes
- Increasing procedural turnaround and the hospital's bottom line

Through our OR Integration Systems, we focus on supporting patients and hospital teams alike, now and into the future, through easy-to-use and scalable technology that supports the operative workflow.



Turning Information into Action

The Harmony iQ 3600 Integration System's focus is clear: deliver critical information – from many video/data, EMR and captured media – when and where it is needed for the OR staff to make timely, accurate patient care decisions, and document the procedure with ease and accuracy.

Scalable to Accommodate up to 32 Connections as Your Needs Grow

Support high-volume rooms, visually complex procedures and demanding hybrid ORs. Connect up to 32 devices and route up to 32 destinations for HD viewing, capturing and recording.

The Harmony iQ 3600 Integration System connects teams within and beyond the OR through robust collaboration and communication tools.




A User-Focused Touch Panel Designed by Nurses, for Nurses

It's critical to ensure that OR teams can use the technology in the OR easily and efficiently. Our control panel was designed by those who use it every day, OR nurses. Operate audio sources, in-light cameras, wall cameras and the Free5® Camera System simply and conveniently. Support efficient and safe procedures, and respond quickly to events that create a need for adjustment to these devices.

More Information. Less Clutter.

Our base unit mounts flush within wall studs for a zero-footprint design, and allows you to connect virtually any imaging or data device from any vendor through a vendor-neutral integration approach. The unit and cabling are kept off the floor, reducing clutter and trip hazards.




Truly Connected for Better Patient Outcomes

Compare and Contrast Video at a Glance

Experience enhanced visualization that allows OR teams to view video from multiple sources side-by-side on one screen.

Rotate images to portrait mode to get a different view of spine cases or other procedures requiring vertical viewing.



Capture Everything. Keep What You Want.

Our **VuCapture® Digital Recording System** enables you to capture surgical procedures in stunning high-definition detail and associate them with patient case* data, then print or save to a USB drive.

Through our **Harmony iQ Perspectives® Image Management System**, easily send captured media to the Electronic Medical Record, and securely review the media from a computer or mobile device. Images and recordings are automatically saved to the correct case to support accurate and thorough case records.*

* Images are not suitable for diagnostic applications.

Happy Patient. Happy Staff.

Through our **SignatureSuite™ OR Integration System** offerings, we can help reduce the patient's anxiety before a procedure by dimming the lights and displaying calming media on monitors, such as a beach scene or mountain landscape as the patient enters the OR. This feature also helps to reduce stress for the OR team and makes set-up run more efficiently.

It's one more tangible way your hospital can demonstrate its focus on the patient. Once the patient is anesthetized, video routing presets exit from SignatureSuite OR Integration System and prepare for the particular type of procedure's visualization requirements.



Connect and Collaborate in Real-Time.

Connect teams within the OR with clinicians and audiences outside the OR, facilitating real-time collaboration and education.

Our **Call and Connect** streaming and conferencing technologies are supported on the Harmony iQ 3600 Integration System, delivering HD audio and video streaming and bi-directional audio and videoconferencing capabilities within the hospital network or beyond.


Ready for 4K? We Are.



4K Planning With Ease and Flexibility

Making your OR 4K compatible can be as easy as upgrading your existing integration equipment.

STERIS' 4K upgrade for the Harmony iQ 2800 Integration System consists of 4K routing hardware that installs with the existing system and the Vividimage® 4K Surgical Displays and Vividimage® W 4K Wall Displays.



Vendor-Neutral OR Integration: Get the Most Out of Your 4K Investment

What good is a 4K integration system that ties your hands to limited imaging devices and capabilities? Our vendor-neutral architecture empowers you to choose the 4K imaging devices that best suit your needs.

Supporting Your Evolving Technology Needs ... at Your Pace

Our scalable design and fiber-optic architecture support more connections and higher resolutions as your needs grow. Through our future-ready design, you can be assured your OR Integration System will meet your needs today and in the future.

Advanced Visualization Without Compromise

We recognize that 4K OR integration is about more than just seeing ultra high-resolution surgical video. With the Harmony iQ 3600 Integration System, you also retain the robust collaboration and workflow tools you expect in a STERIS integrated OR.



Product Specifications



Technology and Support You Can Count On

Skilled STERIS technicians are available by phone and on-site to support real-time identification and resolution of technical issues, as well as perform regular maintenance to keep the system performing at its best. Advanced remote diagnostic support lets STERIS perform many support functions without needing to enter the OR.

We'll quickly diagnose and resolve system issues and help new users learn the ropes, keeping your focus on the patient.

STERIS – A History of OR Innovation

STERIS OR Integration Systems support positive clinical outcomes and efficiency, building on generations of surgical solutions that give OR staff the freedom and confidence to focus on what they do best – care for patients.

For more information, contact your STERIS representative or visit www.steris.com.



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
800-548-4873
www.steris.com