

**APPLICATION**

The Reliance Endoscope Transport System is intended for use as a system consisting of a transport cart, endoscope transport bins and transport bags for transporting ready-to-use and post-procedural flexible endoscopes within a facility.

**DESCRIPTION**

The transport cart contains six shelves for holding six transport bins.

The transport bins are clear/translucent so endoscopes can be easily seen in the transport bags.

The transport bags are supplied in sterile packaging. Each package contains one green transport bag and one red transport bag. The green transport bag is used for ready-to-use endoscopes. The red transport bag is used for post-procedural endoscopes.

**SIZE**

- **Transport Cart (RETCART01)**  
 46-1/2" (1181 mm) – Height  
 27" (686 mm) – Front to Back  
 19-3/4" (502 mm) – Left to Right
- **Transport Bin (RETBIN01)**  
 18 x 24 x 6" (457 x 610 x 152 mm)
- **Transport Bag (RETBAG01)**  
 33 x 21-1/4" (840 x 540 mm)  
**Sterile Packs**  
 5-3/8 x 3" (136 x 76 mm)

**CONSTRUCTION**

The transport cart is made of stainless steel. The transport cart is equipped with medical grade locking/swivel casters.

The transport bins are made of plastic.



Typical only – details may vary

**ACCESSORIES**

- Transport Bin with Lid
- Transport Bags – Case of 200 pre-packaged sterile packs. Each pack contains one green bag and one red bag.

**The base language of this document is ENGLISH. Any translations must be made from the base language document.**

**The Selections Checked Below Apply To This Equipment**

**CAPITAL EQUIPMENT**

- Reliance Endoscope Transport System (RETCART01)

**ACCESSORIES**

- Transport Bin with Lid (RETBIN01)
- Transport Bags (RETBAG01)

Item \_\_\_\_\_

Location(s) \_\_\_\_\_

---

**For Further Information, contact:**



STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 • USA  
440-354-2600 • 800-548-4873  
[www.steris.com](http://www.steris.com)