

Celerity™ 20 HP Biological Indicator

LESS WAITING
LESS WORRYING

MORE CONFIDENCE



Meet OR needs with confidence and on-time with the **FASTEST** BI results.


SATISFYING THE NEED FOR SPEED



Adding the **fastest BI results** helps you deliver sets to the OR on time and with confidence.
NO MORE DELAYS or cancellations to the hospital's busy schedule waiting on BI incubation.

20-MINUTE BI read for **ALL** models and cycles of V-PRO® Low Temperature Sterilizers and the STERRAD® NX® System*, STERRAD® 100NX System® and the STERRAD® 100S® System.

NO MORE WAITING ON RESULTS
NO ONE WAITING ON YOU




Now you can **stay on schedule** and focused on patient care with timely load release.





INNOVATIVE BI DESIGN

Glass-free design helps ensure staff safety



QUICK & EASY

Eliminate the need for additional steps. Simple twist and flick releases the media for incubation.





Speed. Ease of use. Increased productivity.

Reprocess in a V-PRO® Low Temperature Sterilizer and experience accelerated load release. Reducing incubation time to 20 minutes maximizes productivity and minimizes doubt — helping the SPD calmly keep pace with even the busiest OR schedules.

Sterilizer	Sterilizer Cycle	Cycle Time	BI Incubation and Readout	Total Load Release
V-PRO® maX Low Temperature Sterilizer	Non Lumen Cycle	28 minutes	20 minutes	48 minutes
	Flexible Cycle	35 minutes		55 minutes
	Lumen Cycle	55 minutes		75 minutes
V-PRO® 60 Low Temperature Sterilizer	Non Lumen Cycle	28 minutes		48 minutes
	Flexible Cycle	38 minutes		58 minutes
	Lumen Cycle	60 minutes		80 minutes
V-PRO® 1 Plus Low Temperature Sterilizer	Non Lumen Cycle	28 minutes		48 minutes
	Lumen Cycle	55 minutes		75 minutes
V-PRO® 1 Low Temperature Sterilizer	Lumen Cycle	55 minutes		75 minutes
STERRAD® 100 NX	Standard Cycle	47 minutes		67 minutes
	FLEX Cycle	42 minutes		62 minutes
	Express Cycle	24 minutes		44 minutes
	DUO Cycle	60 minutes		80 minutes
STERRAD® NX	Standard Cycle	28 minutes		48 minutes
	Advanced Cycle	38 minutes		58 minutes
STERRAD® 100S	Standard Cycle	55 minutes		75 minutes

The TOTAL package

COMPLETE
sterility assurance solutions
for V-PRO sterilizers.

	PRODUCT NAME	SKU	QUANTITY
	celerity™ 20 HP Biological Indicator	LCB044	25/box
	celerity™ HP Incubator	LCB046	1 each

BECAUSE YOU NEED TO
KNOW FAST

VERIFY™ VH202 Indicator Tape

Immediate confirmation
of sterilant exposure.



VERIFY™ HPU Chemical Indicators

Immediate confirmation of
sterilant penetration.

Vis-U-All™ Low Temperature Sterilization Pouches

Immediate identification of device
and chemical indicator results.



IT'S ABOUT **TIME**

Solutions across
sterilization modalities



No matter which modality, only STERIS offers
biological indicator results in **40 minutes or less**.

Steam Sterilization
40 minutes



Hydrogen Peroxide
Sterilization*
20 minutes



*Current FDA clearance applies to V-PRO® 1 Low Temperature Sterilizer, V-PRO® 1 Plus Low Temperature Sterilizer, V-PRO® maX Low Temperature Sterilizer, V-PRO® 60 Low Temperature Sterilizer, STERRAD® 100NX®, STERRAD® NX® and STERRAD® 100S® Sterilizers.

Visit steris.com/celerity or contact your local
Sales Representative for more information.

 **STERIS®**

STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 • USA
440-354-2600 • 800-548-4873
www.steris.com

STERRAD is a registered trademark of Advanced Sterilization
Products, a division of Johnson and Johnson Company

M9189EN.2017-12, REV.A ©STERIS Corporation.
All Rights Reserved.