



Frontal Sitting Attachment #1802

Frontal Sitting Attachment #1802 mounts neurosurgical headrests to head end of AMSCO 1080 & AMSCO 2080 Surgical Table's back section and head or foot end of AMSCO 3080 & AMSCO 3085 Surgical Tables.

Description

Order Number: BF01

Mounts neurosurgical headrests to head end of Amsco® 1080 & Amsco® 2080 Surgical Table's back section and head or foot end of Amsco® 3080 & Amsco® 3085 Surgical Tables. Ball joints allow full angular positioning of headrest.

Other Attachments Required:

- Order Number: BF03-800, Multi-Poise Headrest or
- Order Number: BF03-600, Horseshoe Headrest
- Order Number: BF18-500, Eye-ENT-Neuro Headrest Adaptor when used with Amsco 3080 and Amsco 3085 Surgical Tables

Table Compatibility:

- Amsco® 1080 Surgical Table
- Amsco® 2080 Surgical Table
- Amsco® 3080 Surgical Table
- Amsco® 3085 Surgical Table
- Amsco® Ortho II Surgical Table
- Amsco® OrthoVision® Orthopedic and Fracture Table

Usage: For support of patient's head from behind when patient is in a seated position for temporal, parietal, and frontal surgical procedures position for temporal, parietal, and frontal surgical procedures.

NOTE: Patient weight limitation is 300 pounds (136 kg).

Literature

Related Products

