



STERIS® 4085 General Surgical Table

The STERIS 4085 General Surgical Table's design gives a surgical team the flexibility, control and convenience they need for today's patients and procedures.

Description

Outstanding versatility and value

The STERIS 4085 General Surgical Table's design gives a surgical team the flexibility, control and convenience they need for today's patients and procedures.

Features for a wide array of procedures:

- Accommodates virtually all patients with generous weight and height ranges
- Provides outstanding access for C-arm, permitting high quality images crucial to cardiothoracic, orthopedic and neurosurgical procedures
- Contours to virtually all patient postures with four section tabletop
- Compensates for uneven OR floors with self-leveling floor locks
- Radiolucent 18" sliding top for total body imaging (43" upper and 43" lower body C-arm access)
- Provides unrestricted perineal access with removable leg section
- Optional Featherweight Leg Section complies with [AORN Safe Handling Principles](#)
- Enhances lateral procedures with a powered, radiolucent kidney elevator
- Easy to use, easy to read hand control with digital data window
- Prevents sectional and table base collisions with the Auto Limit Sensor™
- Compatible with a wide array of previously purchased accessories

Description

Outstanding versatility and value

The STERIS 4085 General Surgical Table's design gives a surgical team the flexibility, control and convenience they need for today's patients and procedures.


Features for a wide array of procedures:

- Accommodates virtually all patients with generous weight and height ranges
- Provides outstanding access for C-arm, permitting high quality images crucial to cardiothoracic, orthopedic and neurosurgical procedures
- Contours to virtually all patient postures with four section tabletop
- Compensates for uneven OR floors with self-leveling floor locks
- Radiolucent 18" sliding top for total body imaging (43" upper and 43" lower body C-arm access)
- Provides unrestricted perineal access with removable leg section
- Optional Featherweight Leg Section complies with [AORN Safe Handling Principles](#)
- Enhances lateral procedures with a powered, radiolucent kidney elevator
- Easy to use, easy to read hand control with digital data window
- Prevents sectional and table base collisions with the Auto Limit Sensor™
- Compatible with a wide array of previously purchased accessories

Positioning



Femoro-Popliteal



4085 : Cardiovascular

Cardiovascular



Cystoscopy/Urology/GYN



Endourology



Ophthalmic/ENT



Bariatric/Split Leg



Kidney & Thoracic



Orthopedic



Shoulder Chair



Total Hip



Prone Neurosurgery



Spinal Neurosurgery



Lateral Neurosurgery



Kneeling Neurosurgery



Nisson Fundoplication



Gall Bladder and Thoracic



Endovascular



Gynecology

Specifications

SPECIFICATION	DESCRIPTION
Table Length	81" (2057 mm)
Table Width	20" (508 mm) (28° with extenders)
Table Height Range	26" to 45" (660 to 1143 mm)
Table Slide Range	9" (227 mm) to head, 9" to foot (18" total)
Patient Weight Capacity	1,100 lbs. (500 kg) patient support, including raise (centered on the column) 1,000 lbs (454 kg) full table articulation (centered column) 600 lbs (272 kg) full table articulation, including slide
Table Shipping Weight	560 lbs. (254 kg)
Trendelenburg / Reverse	30° / 30°
Lateral Tilt (left/right)	20° / 20°
Head Section	+90° to -90°
Back Section	+80° to -40°
Leg Section	0° to -105° (removable)
Flex/Reflex	140° / 100°
Manual Override	Yes
Perneal Cut-out	Yes
Kidney Elevator	4" powered radiolucent

Literature

 [Adobe® Acrobat® Reader required](#)

BROCHURE

Document #





Document Title

 [M3309EN](#)





STERIS 4085 SURGICAL TABLE: A SMART INVESTMENT IN SURGICAL VERSATILITY

 [M4038FC](#)

SECURECARE TABLE REVIVE BROCHURE-FRENCH CANADIAN VERSION

-  [M4038EN](#) SECURECARE TABLE REVIVE BROCHURE
-  [M2699EN](#) POWER-LIFT AND BARIATRIC POWER-LIFT STIRRUPS. INFINITE LEG POSITIONING. PRECISELY.
-  [M3483EN](#) ADDRESSING A PRESSING NEED. STERIS HLT SURGICAL TABLE PAD
-  [M4099EN](#) SECURECARE TABLE CARE PACKAGE SELL SHEET

SELL SHEET

Document #	Document Title
 M4099FC	SECURECARE TABLE CARE PACKAGE SELL SHEET
 M4060EN	SURE LOC XPS SELL SHEET
 M4034EN	NEW SURGICAL TABLE ACCESSORIES
 M4061EN	E-Z LIFT BEACH CHAIR SELL SHEET

Related
Products

Surgical Table Accessories



Anesthesia Armboard

- Lightweight, carbon fiber design
- Distal end trigger for full horizontal rotation
- Universal compatibility to North American standard siderails



Power-lift Stirrups

- Ergonomic, pneumatic lift-assist for safe patient handling
- Single handle action for full articulation and range of motion



Shoulder Chair

- Used for open or arthroscopic procedures
- Features a unique, multi-axis head positioning capability
- Includes head restraint device and integral lateral brace



Pad Positioners

- Used in supplement with tabletop pads for promotion of proper patient positioning
- Provide stability and pressure management

Surgical Tables



STERIS 5085 SRT Surgical Table

- Slide, 360° Rotate and Transport capabilities
- Largest clinically-relevant image area
- Generous patient height and weight ranges in a sliding top design for improved imaging



STERIS 5085 Surgical Table

- Low height for performance suites requiring versatility
- Provides unrestricted perineal access with a removable leg section

Operating Room Equipment



Surgical Lights

- Surgical, exam, and labor and delivery lighting systems
- Deliver pure white light, and enhancement of deep saturated reds and colors



Surgical Booms

- Sleek, compact design to allow nurse line of sight to the sterile field
- Straight drop units and single or dual arm units



Surgical Displays

- Surgical field and wall-mounted displays
- Deliver HD and 4K visualization within the surgical field



OR Integration

- Designed to improve workflow, procedure guidance and peer collaboration
- Deliver scalable solutions to best fit the diverse needs of the OR

Service
Manual

Service Manual



[STERIS 4085 General Surgical Table
Maintenance Manual](#)

Replacement
Parts

Genuine STERIS Parts for Your STERIS Surgical Tables



[View Service Parts](#)



By choosing genuine STERIS OEM replacement parts, you can be assured the part is chosen based on the highest quality standards to ensure your equipment is the best. STERIS's OEM replacement parts promote compliance with regulatory and organizational standards on equipment maintenance and repair.