



STERIS® 5085 General Surgical Table

The STERIS 5085 General Surgical Table is a mobile, electrohydraulically operated surgical table that delivers with industry-leading innovations, specifically designed for the OR of today and tomorrow.

Description

The STERIS 5085 General Surgical Table quickly and easily adapts to virtually every surgical specialty. Higher weight limit, exceptionally low height and a large imaging area means no more inflexible, dedicated specialty tables.

Quick, one-touch patient posturing is enabled by the 5085 Surgical Table's intuitive hand control's pre-programmed positions for rotation and beach chair for easy set up. The surgical table's easy-to-read LCD window telegraphs table positioning data to caregivers, even in a darkened room.

The STERIS 5085 Surgical Table:

- Accommodates virtually all patients with generous weight and height ranges
- Provides unrestricted perineal access with removable leg section
- Optional Featherweight Leg Section complies with AORN Safe Handling principles
- Enhances lateral procedures with a powered, radiolucent kidney elevator
- Provides outstanding access for C-arm, permitting clear high-quality images crucial to MIS, cardiothoracic, orthopedic and neurosurgical procedures
- Facilitates quick and safe patient positioning with an easy-to-use, easy-to-read hand control
- Contours to virtually all patient postures with four section tabletop
- Compensates for uneven OR floors with self-leveling floor locks
- Prevents sectional and table base collisions with the Auto Limit Sensor
- Compatible with a wide array of previously purchased accessories

Flexible Design

The 5085 Surgical Table's strong, yet flexible design provides an exceptionally low height for MIS, robotic, and bariatric techniques without sacrificing the height needed for orthopedic and cardiac procedures.

Clear Imaging

The surgical table's large imaging area permits clear, high quality images crucial for today's MIS, cardiothoracic, orthopedic and neurological procedures, while the narrow column and generous longitudinal slide allows 100% access for the C-arm and surgical team without patient reorientation.

Intuitive Controls

The 5085 Surgical Table's user-friendly, data-rich hand control offers:

- **An easy-to-read LCD window** – provides information on table positioning, service requirements, and potential interference
- **Pre-programmed, backlit buttons** – quick, one-touch patient posturing
- **4-section tabletop** – articulates to virtually every patient posture
- **LEDs** – highlight key buttons such as floor lock, normal and reverse patient orientation
- **Touch pads** – include menu advancement, programmed positions, and more

Featherweight Leg Section

The modular, removable leg section's carbon fiber construction makes it strong, lightweight and ideal for imaging. The surgical table's featherweight leg section allows for flexibility in movement and function to experience standard articulation capabilities.

The STERIS 5085 General Surgical Table quickly and easily adapts to virtually every surgical specialty. Higher weight limit, exceptionally low height and a large imaging area means no more inflexible, dedicated specialty tables.

Quick, one-touch patient posturing is enabled by the 5085 Surgical Table's intuitive hand control's pre-programmed positions for rotation and beach chair for easy set up. The surgical table's easy-to-read LCD window telegraphs table positioning data to caregivers, even in a darkened room.

The STERIS 5085 Surgical Table:

- Accommodates virtually all patients with generous weight and height ranges
- Provides unrestricted perineal access with removable leg section
- Optional Featherweight Leg Section complies with AORN Safe Handling principles
- Enhances lateral procedures with a powered, radiolucent kidney elevator
- Provides outstanding access for C-arm, permitting clear high-quality images crucial to MIS, cardiothoracic, orthopedic and neurosurgical procedures
- Facilitates quick and safe patient positioning with an easy-to-use, easy-to-read hand control
- Contours to virtually all patient postures with four section tabletop
- Compensates for uneven OR floors with self-leveling floor locks
- Prevents sectional and table base collisions with the Auto Limit Sensor
- Compatible with a wide array of previously purchased accessories

Flexible Design

The 5085 Surgical Table's strong, yet flexible design provides an exceptionally low height for MIS, robotic, and bariatric techniques without sacrificing the height needed for orthopedic and cardiac procedures.

Clear Imaging

The surgical table's large imaging area permits clear, high quality images crucial for today's MIS, cardiothoracic, orthopedic and neurological procedures, while the narrow column and generous longitudinal slide allows 100% access for the C-arm and surgical team without patient reorientation.

Intuitive Controls


The 5085 Surgical Table's user-friendly, data-rich hand control offers:

- **An easy-to-read LCD window** – provides information on table positioning, service requirements, and potential interference
- **Pre-programmed, backlit buttons** – quick, one-touch patient posturing
- **4-section tabletop** – articulates to virtually every patient posture
- **LEDs** – highlight key buttons such as floor lock, normal and reverse patient orientation
- **Touch pads** – include menu advancement, programmed positions, and more


Featherweight Leg Section

The modular, removable leg section's carbon fiber construction makes it strong, lightweight and ideal for imaging. The surgical table's featherweight leg section allows for flexibility in movement and function to experience standard articulation capabilities.


Positioning




Femoro-Popliteal




Cardiovascular



Cystoscopy/Urology/GYN




Sitting Neurosurgery




Ophthalmic/ENT




Bariatric/Split Leg



Kidney & Thoracic



Orthopedic



Shoulder Chair

Other Positions

Total Hip
Prone Neurosurgery
Spinal Neurosurgery
Endovascular

Lateral Neurosurgery
Kneeling Neurosurgery
Nisson Fundoplication





Lithotomy
Total Hip
Gall Bladder & Thoracic

Specifications






Literature

 [Adobe® Acrobat® Reader required](#)



BROCHURE

Document #	Document Title
 M3302EN	STERIS 5085 AND 5085 SRT SURGICAL TABLES: A NEW BENCHMARK FOR SURGICAL VERSATILITY
 M2699EN	POWER-LIFT AND BARIATRIC POWER-LIFT STIRRUPS. INFINITE LEG POSITIONING. PRECISELY.
 M3483EN	ADDRESSING A PRESSING NEED. STERIS HLT SURGICAL TABLE PAD
 M4099EN	SECURECARE TABLE CARE PACKAGE SELL SHEET

SELL SHEET

Document #	Document Title
 M9750EN	SHORT FEATHERWEIGHT LEG SECTION SELL SHEET
 M4099FC	SECURECARE TABLE CARE PACKAGE SELL SHEET
 M3355EN	THE ONLY STERIS-APPROVED DISINFECTANTS/CLEANERS FOR YOUR SURGICAL TABLE PADS
 M4060EN	SURE LOC XPS SELL SHEET
 M4034EN	NEW SURGICAL TABLE ACCESSORIES

TECH DATA

Document #	Document Title
 10270850	STERIS 5085 GENERAL SURGICAL TABLE
 SD901	STERIS 5085 GENERAL SURGICAL TABLE

Related
Products

Replacement
Parts

Service
Manual