



SYSTEM 1[®] endo Liquid Chemical Sterilant Processing System

- For safe and effective reprocessing of cleaned, reusable, heat-sensitive, **semi-critical** medical devices and their accessories
- Rapid 18-minute cycle
- Eliminates all microbial life, including bacterial endospores, by utilizing S40[®] Sterilant Concentrate

Description

The SYSTEM 1 endo Liquid Chemical Sterilant Processing System is a liquid chemical sterilant processor designed to safely and effectively reprocess heat-sensitive, semi-critical medical devices.

How the SYSTEM 1 endo Liquid Chemical Sterilant Processing System Works

The SYSTEM 1 endo Processor ensures that S40 Sterilant use dilution reaches all parts of the endoscope, including port surfaces commonly shielded from germicide contact in Automated Endoscope Reprocessors (AERs) with tight-fitting adapters. Quick Connects provide flow of chemistry throughout the endoscope lumens and around adapter ports.

Why the SYSTEM 1 endo Liquid Chemical Sterilant Processing System?

The SYSTEM 1 endo Processor is validated for efficacy and material compatibility for over 2,000 devices including colonoscopes, gastroscopes, bronchoscopes, enteroscopes, duodenoscopes and more. The validated process of the SYSTEM 1 endo Processor goes beyond high-level disinfection to eliminate all microbial life, including bacterial spores. This reduces the risk of microbial transmission by processed medical devices.

Device Compatibility

Device Compatibility Matrix

Confirm that your medical devices and accessories are compatible for processing with our Low Temperature Systems

[search for my devices](#)

Literature

BROCHURE

Document #	Document Title
	SYSTEM 1 ENDO BROCHURE

TECH DATA

Document #	Document Title
	SYSTEM 1 ENDO LIQUID CHEMICAL STERILANT PROCESSING SYSTEM

TECHNICAL MONOGRAPH

Document #	Document Title
	TECHNICAL DATA MONOGRAPH FOR SYSTEM 1 ENDO LIQUID CHEMICAL STERILANT PROCESSING SYSTEM

WALL CHARTS

Document #	Document Title
	PROCESSOR WALL CHART FOR SYSTEM 1 ENDO LIQUID CHEMICAL STERILANT PROCESSING SYSTEM

Related Products

Sterilant



S40® Sterilant Concentrate

- Designed for use in STERIS automated Liquid Chemical Sterilant Processing Systems

Sterility Assurance



VERIFY® Chemical Indicator for S40® Sterilant

- Verifies that the minimum recommended concentration of S40 Sterilant Concentrate use dilution has been met.
- Designed for use with STERIS automated Liquid Chemical Sterilant Processing Systems.



VERIFY® Spore Test Strip for S40® Sterilant

- Designed for use with S40 Sterilant
- Tests sporidical activity of the S40 Sterilant

Printer Paper and Ribbon



Printer Paper

- Box of three rolls



Printer Ribbon

- Box of two cartridges